

### Scientific Products

The Clements Coombs Test Centrifuge is specially designed to prepare washed red cells for the negative Coombs Test.



proven through performance

- **Repeatability** Centrifuge operates for exactly the time set. Assists in eliminating false readings.
- **Secured Tubes** Removable head allows loading, conveying, and decanting to take place **without** having to remove the tubes from the head
- **Stainless Steel Bowl** Bowl unaffected by saline or blood spillage.
- **Safety** Lid interlocked power cut off switch

The Coombs Test centrifuge is also suitable for blood grouping, typing and cross matching, spinning down samples of whole blood, urines and other body fluids.

### **Ordering**

**CEN 96225** Coombs Test Centrifuge with Digital Timer complete with CEN 96242 Fixed Angle 12 Place Head

### **Technical Specifications**

<b>Bowl</b>	Stainless steel
<b>Timer</b>	4 digit, 0-99 mins 59 secs
<b>End of Cycle</b>	Audible alert
<b>Indicators</b>	Min, Sec, Start, Stop 4 Digit LCD time display Run LED
<b>Power</b>	220V/240V - 50 Hz.
<b>Speed</b>	2800 rpm
<b>R. C. F.</b>	700 ( with Cat. No. CEN 96242 Head)
<b>Safety</b>	Lid interlocked power cut off switch
<b>Unit Dimensions</b>	Height 245 mm, Width 350 mm, Depth 370 mm
<b>Shipping Dimensions</b>	Height 350 mm, Width 400 mm, Depth 450 mm
<b>Shipping Weight</b>	18 kg

\*Specifications subject to change without notice.

**The Clements range of medical and scientific products is renowned worldwide for quality, reliability and durability**

Distributed By:

Proudly designed and manufactured in Australia by:

**Clements Medical Equipment**  
5 Cal Close Somersby NSW 2250  
Australia  
Tel +61 2 4340 2484  
Fax +61 2 4340 5121  
Email: [sales@clements.net.au](mailto:sales@clements.net.au)  
Website: [www.clements.net.au](http://www.clements.net.au)

