

# Orbital 325 Benchtop Centrifuge

The Clements Orbital 325 is designed for clinical and small laboratory centrifugation where uncomplicated analogue speed and digital timing control is required. The Orbital 325 is available with either a six place rotor for sealed tube carriers or an eight place rotor for open tube carriers.

proven through performance



- **Functional Controls** For easy and precise operation
- **Accurate** Digital timer for precise timed operations
- **Electric Lid Lock** Lid lock prevents accidental opening of the centrifuge while running
- **Economical** Provides basic centrifuge functions for standard requirements
- **Aerosol Protection** Aerosol sealed tube carriers (CEN 91517)

<b>ORDERING</b>	
<b>CEN 91516</b>	Orbital 325 Centrifuge with Integrated Lid Lock, Digital Timer and 8 Place Rotor for Open Tube Carriers (Carriers Supplied)
<b>CEN 91517</b>	Orbital 325 Centrifuge with Integrated Lid Lock, Digital Timer and 6 Place Rotor for Sealed Tube Carriers (Carriers Supplied)
<b>TECHNICAL SPECIFICATIONS</b>	
<b>Capacity</b>	8 x 15 ml Open Tube Carriers (CEN 91516) 6 x 15 ml Sealed Tube Carriers (CEN 91517)
<b>Speed</b>	Variable speed control from 1000 to 3200 RPM
<b>RCF (max)</b>	1,700 g
<b>Digital Timer</b>	Countdown 0-99 mins 59 secs
<b>Speed Control</b>	Analog Dial
<b>Indicators</b>	Min, Sec, Start, Stop 4 Digit LCD time display Run LED
<b>Safety</b>	Integrated Electronic Lid Lock
<b>Power</b>	240 V 50 Hz 0.5A
<b>Shipping Weight</b>	16 kg
<b>Unit Dimensions</b>	Height 245 mm, Width 350 mm, Depth 370 mm
<b>Shipping Dimensions</b>	Height 350 mm, Width 400 mm, Depth 450 mm

All specifications are subject to change without notice.  
Please check the exact specifications when ordering.

**The Clements range of medical and scientific products is renowned worldwide for quality, reliability and durability.**

Distributed by:

Proudly designed and manufactured in Australia by:

**Clements Medical Equipment Pty Ltd**

5 Cal Close Somersby NSW 2250

Australia

Tel: +61 2 4340 2484

Fax: +61 2 4340 5121

Email: [sales@clements.net.au](mailto:sales@clements.net.au)

Web site: [www.clements.net.au](http://www.clements.net.au)

