

| Ordering | | | | |
|------------------------------|-----------|-----------|-----------|-----------|
| Speed | 6 RPM | 12 RPM | 20 RPM | 30 RPM |
| Rotor Capacity | | | | |
| 18 x 15 ml | CEN 96151 | CEN 96152 | CEN 96153 | CEN 96154 |
| 9 x 15 ml + 9 x 50 ml | CEN 96161 | CEN 96162 | CEN 96163 | CEN 96164 |
| 18 x 50 ml | CEN 96171 | CEN 96172 | CEN 96173 | CEN 96174 |

Other rotor clip and motor speed variations are available on request Please contact Clements for further information.

| Technical Specifications | |
|----------------------------|--|
| Switch Controls | On / Off |
| Materials: | |
| Disc | Satin stainless steel |
| Clips | Phosphor bronze, chrome plated, double riveted |
| Drive | Friction drive centre shaft |
| Overload Protection | Safety drive clutch |
| Stability | Load insensitive gearbox design |
| Power | 240 V 50 Hz 0.2 A |
| Dimensions (mm) | 270 H • 260 W • 248 D |
| Weight | 3.2 kg |

The Clements range of medical and scientific products is renowned worldwide for quality, reliability and durability

*Specifications subject to change without notice.

| |
|-----------------|
| Distributed By: |
|-----------------|

| | |
|--|--|
| <p>Proudly designed and manufactured in Australia by: Clements Medical Equipment Pty Ltd 5 Cal Close, Somersby NSW 2250 Australia</p> <p>Tel +61 2 4340 2484 Fax +61 2 4340 5121 Email: sales@clements.net.au Web site: www.clements.net.au</p> |  <p>AUSTRALIAN MADE[®]</p> |
|--|--|