

SUCTION CANISTERS

DISPOSABLE LINER SYSTEM

Clements offers the FLOVAC range of suction collection jars / canisters with disposable liners with the following features.

- Available in 1L, 2L and 3L sizes.
- Inbuilt Overflow Protection (non-mechanical)
Hydrophobic filter with anti-microbial protection (bacteria and viruses), volatile particulate and electro-surgery smoke protection.
- Secure handling for improved infection control
- Adaptable quick fit mounting brackets



CLEMENTS

www.clements.net.au

DISPOSABLE LINERS AND ACCESSORIES



DISPOSABLE LINERS AND ACCESSORIES

Group	Description	Catalogue Number	①	②	③	④	⑤	⑥
Liners	Disposable FLOVAC 1L Liner (Pack 50)	SUC81030121				50		
	Disposable FLOVAC 2L Liner (Pack 50)	SUC81030122				50		
	Disposable FLOVAC 3L Liner (Pack 50)	SUC84500123				50		
Support Jars	Jar, Support for FLOVAC 1L Disposable Liner	SUC81030005			1			
	Jar, Support for FLOVAC 2L Disposable Liner	SUC81030132			1			
	Jar, Support for FLOVAC 3L Disposable Liner	SUC84500013			1			
Mounts	Mounting Block for Ring Bracket, 25x5	SUC84500031	1					
	Ring Bracket for FLOVAC Support Jar	SUC84500032		1				
	Nipple, Yellow for FLOVAC Disposable Liner	SUC84500033						2

"VACUUM" Connector	Female conical to suit SUC84500033
"PATIENT" Connector	Ø 14.0- 15.5 mm (Ø 8.0- 9.2 mm with supplied elbow)
"TANDEM" Connector	Ø 8.0- 9.2 mm
Maximum Vacuum	-95 kPa / 5 minutes
Recommended Patient Tubing Size	Ø 6.0 x 9.0 mm, L = 2.5 m max
Support Jar Cleaning	Autoclave 121°C for 15 min
Support Jar Material	Polycarbonate
Liner Bag Material	LDPE (low density polyethylene)
Liner Lid Material	HDPE (high density polyethylene)

All specifications are subject to change without notice. Please check the exact specification when ordering.

ICU Medical Australia
Unit U, 10- 16 South St
Rydalmere NSW 2116, Australia
Telephone +61 2 9466 5300

Rev 5 © Copyright ICU Medical Australia 2023

CLEMENTS